


Anti-Tom70 Rabbit Monoclonal Antibody For Wb, Elisa - 094826

Item Number: CM0002071



Introduction

High-specificity rabbit anti-TOM70 monoclonal antibody validated for WB and ELISA.

Recognizes human, mouse, and rat TOM70, a mitochondrial import receptor and antiviral response mediator. Ideal for mitochondrial research and viral infection studies.

[Learn More](#)

Parameter	Value
Product Name	TOM70 Rabbit mAb
Remarks/Alias	Tom70; TOMM70A; TOM70
Species	Human
Gene ID (Human)	9868
Gene ID	9868
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 245-608 of human TOM70 (NP_055635.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O94826
Protein Weight	67kDa
Shipping	Ice bag