

Anti-Tom70 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O94826

Item Number: CM0028548



Introduction

Rabbit polyclonal antibody targeting human TOM70 (mitochondrial import receptor subunit TOM70). Validated for WB, IF/ICC, and ELISA; cross-reacts with mouse and rat. Ideal for studying mitochondrial protein import and antiviral signaling.

[Learn More](#)

Parameter	Value
Product Name	TOM70 Rabbit pAb
Remarks/Alias	Tom70; TOMM70A; TOM70
Species	Human
GeneID (Human)	9868
GeneID	9868
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 245-608 of human TOM70 (NP_055635.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O94826
Protein Weight	67 kDa
Shipping	Ice bag