

Anti-Tomm22 Rabbit Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q9Ns69

Item Number: CM0011027



Introduction

High specificity rabbit monoclonal antibody targeting human TOMM22, validated for Western blot, IHC-P, IP, and ELISA. Recognizes mitochondrial import receptor subunit TOM22 homolog across human, mouse, and rat. Ideal for mitochondrial research.

[Learn More](#)

Parameter	Value
Product Name	TOMM22 Rabbit mAb
Remarks/Alias	1C9-2; TOM22; MST065; MSTP065; TOMM22
Species	Human
GeneID (Human)	56993
GeneID	56993
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NS69
Protein Weight	16kDa
Shipping	Ice bag