

Anti-Tomm34 Rabbit Monoclonal Antibody - Q15785

Item Number: CM0010621

Introduction

High-quality rabbit monoclonal antibody targeting TOMM34, a mitochondrial import receptor. Validated for Western blot, IHC-P, and ELISA with human and rat cross-reactivity. Ideal for studying mitochondrial protein import and precursor chaperone functions.

[Learn More](#)

Parameter	Value
Product Name	TOMM34 Rabbit mAb
Remarks/Alias	TOM34; URCC3; HTOM34P
Species	Human
Gene ID (Human)	10953
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-309 of human TOMM34 (NP_006800.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q15785
Protein Weight	35kDa
Shipping	Ice bag