

## Anti-Tomm40 Monoclonal Antibody For Wb - O96008

Item Number: CM0005356



### Introduction

Anti-TOMM40 rabbit monoclonal antibody, validated for WB, IP, IHC-P, IF/ICC, ELISA. Recognizes human, mouse, rat TOMM40. Use for mitochondrial import and mitophagy research.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | TOMM40 Rabbit mAb                                                                                                   |
| Remarks/Alias    | TOM40; PEREC1; C19orf1; PER-EC1; D19S1177E; TOMM40                                                                  |
| Species          | Human                                                                                                               |
| GeneID (human)   | 10452                                                                                                               |
| GeneID           | 10452                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-90 of human TOMM40 (NP_001122388.1) |
| Source           | Rabbit                                                                                                              |
| Category         | Monoclonal Antibodies                                                                                               |
| Application      | WB, IP, IHC-P, IF/ICC, ELISA                                                                                        |
| Cross Reactivity | Human, Mouse, Rat                                                                                                   |
| SWISS            | O96008                                                                                                              |
| Protein Weight   | 38kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |