

Anti-Tomm40 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O96008

Item Number: CM0026621



Introduction

High-quality rabbit polyclonal antibody against TOMM40, a key mitochondrial translocase component. Suitable for WB, IHC-P, IF/ICC, and ELISA applications; cross-reactive with human, mouse, and rat samples. Ideal for mitochondrial import and mitophagy research.

[Learn More](#)

Parameter	Value
Product Name	TOMM40 Rabbit pAb
Remarks/Alias	TOM40; PEREC1; C19orf1; PER-EC1; D19S1177E; TOMM40
Species	Human
GeneID (Human)	10452
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-90 of human TOMM40 (NP_001122389.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O96008
Protein Weight	38kDa
Shipping	Ice bag