

Anti-Tomm40L Rabbit Polyclonal Antibody For Wb, Elisa - Q969M1

Item Number: CM0020995

Introduction

High-quality anti-TOMM40L rabbit polyclonal antibody validated for Western blot (WB) and ELISA applications. Detects human TOMM40L, cross-reacts with mouse. Ideal for mitochondrial protein import research.

[Learn More](#)

Parameter	Value
Product Name	TOMM40L Rabbit pAb
Remarks/Alias	TOMM40B; TOMM40L
Species	Human
GeneID (Human)	84134
GeneID	84134
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human TOMM40L (NP_115550.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q969M1
Protein Weight	33kDa
Shipping	Ice bag