

Anti-Tomm5 Rabbit Polyclonal Antibody For If/Icc And Elisa - Q8N4H5

Item Number: CM0022608



Introduction

Rabbit polyclonal TOMM5 antibody validated for immunofluorescence (IF/ICC) and ELISA. Detects human and mouse TOMM5, a TOM complex subunit involved in mitochondrial protein import and mitophagy signaling.

[Learn More](#)

Parameter	Value
Product Name	TOMM5 Rabbit pAb
Remarks/Alias	Tom5; C9orf105; bA613M10.3; TOMM5
Species	Human
Gene ID (Human)	401505
Immunogen	Synthetic peptide corresponding to amino acids 1-51 of human TOMM5 (NP_001001790.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8N4H5
Protein Weight	6kDa
Shipping	Ice bag