

Anti-Top1Mt Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q969P6

Item Number: CM0024548

Introduction

High-quality rabbit polyclonal antibody against human mitochondrial DNA topoisomerase I (TOP1MT), validated for WB, IHC-P, and ELISA applications. Reacts with human and rat TOP1MT. Suitable for mitochondrial biology research.

[Learn More](#)

Parameter	Value
Product Name	TOP1MT Rabbit pAb
Remarks/Alias	TOP1MT
Species	Human
GeneID (Human)	116447
GeneID	116447
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 470-550 of human TOP1MT (NP_443195.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q969P6
Protein Weight	70 kDa
Shipping	Ice bag