


Anti-Top2A (Dna Topoisomerase Ii Alpha) Rabbit Polyclonal Antibody For Wb And Elisa - P11388

Item Number: CM0020928



Introduction

Rabbit polyclonal antibody against human TOP2A, validated for Western blot and ELISA. Targets DNA topoisomerase II alpha, a key enzyme in DNA topology and chromosome segregation. Ideal for cell cycle and cancer research.

[Learn More](#)

Parameter	Value
Product Name	DNA topoisomerase II alpha (TOP2A) Rabbit pAb
Remarks / Alias	TOP2; TP2A; TOPIIA; TOP2alpha; DNA topoisomerase II alpha (TOP2A)
Species	Human
Gene ID	7153
Immunogen	Synthetic peptide (aa 1100-1200 of human TOP2A, NP_001058.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P11388
Protein Weight	174 kDa
Shipping	Ice bag