


Anti-Topo II α Polyclonal Antibody For Wb, Ihc-P, Elisa - P11388

Item Number: CM0021945



Introduction

High-quality rabbit polyclonal antibody targeting human DNA topoisomerase 2-alpha (Topo II α). Suitable for Western blot, IHC-P, and ELISA. Validated for human samples. Ideal for cancer research and cell cycle studies.

[Learn More](#)

Parameter	Value
Product Name	Topo II α Rabbit pAb
Remarks/Alias	TOP2A; TOP2; TP2A; DNA topoisomerase 2-alpha; Topo II α
Species	Human
GeneID (human)	7153
GeneID	7153
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1431-1531 of human Topo II α (NP_001058.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P11388
Protein Weight	174kDa/177kDa/178kDa/182kDa
Shipping	Ice bag