

Anti-Tor4A Monoclonal Antibody For Wb, Elisa - Q9Nxh8

Item Number: CM0002725

Introduction

Rabbit monoclonal antibody against human TOR4A (Torsin-4A) for Western blot and ELISA applications. Cross-reacts with mouse. Immunogen: recombinant fusion protein spanning amino acids 1-130 of human TOR4A (Q9NXH8).

[Learn More](#)

Parameter	Value
Product Name	TOR4A Rabbit mAb
Remarks/Alias	C9orf167; TOR4A
Species	Human
Gene ID (Human)	54863
Gene ID	54863
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human TOR4A (Q9NXH8)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NXH8
Protein Weight	47kDa
Shipping	Ice bag