

Anti-Torc2/Crtc2 Rabbit Monoclonal Antibody For Wb And Elisa - Q53Et0

Item Number: CM0002143



Introduction

Rabbit monoclonal antibody targeting human TORC2 (CRTC2), validated for Western blot (WB) and ELISA; cross-reacts with rat. Useful for studying CREB-mediated transcription, gluconeogenesis, and mitochondrial biogenesis.

[Learn More](#)

Parameter	Value
Product Name	TORC2 Rabbit mAb
Remarks/Alias	TORC2; TORC-2
Species	Human
Gene ID (Human)	200186
Gene ID	200186
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 594-693 of human TORC2 (Q53ET0).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q53ET0
Protein Weight	73kDa
Shipping	Ice bag