

Anti-Torsina/Tor1A Polyclonal Antibody For Wb, Ihc-P, Elisa - O14656

Item Number: CM0019693



Introduction

TorsinA/TOR1A Rabbit pAb for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Targets the dystonia-related chaperone Torsin-1A involved in protein folding and synaptic vesicle recycling. Ideal for neurobiology research.

[Learn More](#)

Parameter	Value
Product Name	TorsinA/TOR1A Rabbit pAb
Remarks/Alias	DQ2; AMC5; DYT1; TorsinA/TOR1A
Species	Human
GeneID (Human)	1861
GeneID	1861
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 63-332 of human TorsinA/TOR1A (NP_000104.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14656
Protein Weight	38kDa
Shipping	Ice bag