

Anti-Torsina/Tor1A Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - O14656

Item Number: CM0011751



Introduction

Rabbit monoclonal antibody against human TorsinA/TOR1A (O14656), validated for WB, IF/ICC, and ELISA in human, mouse, and rat. This reagent facilitates investigation of protein folding, synaptic vesicle recycling, and nuclear envelope dynamics in dystonia and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	TorsinA/TOR1A Rabbit mAb
Remarks/Alias	DQ2; AMC5; DYT1; TorsinA/TOR1A
Species	Human
GeneID (Human)	1861
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 233-332 of human TorsinA/TOR1A (O14656)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14656
Protein Weight	38kDa
Shipping	Ice bag