

Anti-Tox Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O94900

Item Number: CM0029675

Introduction

TOX rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, rat. Targets thymocyte selection-associated high mobility group box protein TOX, a transcriptional regulator critical for immune cell development and neural differentiation.

[Learn More](#)

Parameter	Value
Product Name	TOX Rabbit pAb
Remarks/Alias	TOX1; TOX
Species	Human
GeneID (Human)	9760
GeneID	9760
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-165 of human TOX (NP_055544.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O94900
Protein Weight	58kDa
Shipping	Ice bag