

Anti-Tp53Aip1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Hcn2

Item Number: CM0029366

Introduction

TP53AIP1 Rabbit Polyclonal Antibody validated for WB and ELISA, with cross-reactivity to mouse and rat. Target molecular weight 13kDa. Suitable for p53-dependent apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	TP53AIP1 Rabbit pAb
Remarks/Alias	P53AIP1; TP53AIP1
Species	Human
Gene ID (Human)	63970
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9HCN2
Protein Weight	13kDa
Shipping	Ice bag