

Anti-Tp53Bp2 Rabbit Polyclonal Antibody For Ihc-P And Elisa - Q13625

Item Number: CM0029711



Introduction

Rabbit polyclonal antibody targeting TP53BP2 (ASPP2), validated for IHC-P and ELISA. Reacts with human, mouse, and rat. Used in apoptosis and cell cycle research. Specifications and technical support available.

[Learn More](#)

Parameter	Value
Product Name	TP53BP2 Rabbit pAb
Remarks/Alias	BBP, 53BP2, ASPP2, P53BP2, PPP1R13A, TP53BP2
Species	Human
Gene ID (Human)	7159
Gene ID	7159
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 360-600 of human TP53BP2 (NP_001026855.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13625
Protein Weight	126kDa
Shipping	Ice bag