

Anti-Tp53Bp2 Rabbit Polyclonal Antibody For Wb, Elisa - Q13625

Item Number: CM0023463



Introduction

Rabbit polyclonal antibody targeting TP53BP2 (ASPP2/53BP2), validated for WB and ELISA in human samples. Regulator of apoptosis and cell growth. Recombinant immunogen, 126kDa. Ideal for apoptosis and cancer research, targeting a key p53-binding protein involved in cell cycle arrest and tumor suppression.

[Learn More](#)

Parameter	Value
Product Name	TP53BP2 Rabbit pAb
Remarks/Alias	BBP; 53BP2; ASPP2; P53BP2; PPP1R13A; TP53BP2
Species	Human
GeneID (Human)	7159
GeneID	7159
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 200-330 of TP53BP2 (NP_001026855.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q13625
Protein Weight	126kDa
Shipping	Ice bag