

Anti-Tp53I13 Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Nbr0

Item Number: CM0020782

Introduction

Rabbit polyclonal antibody against human TP53I13, validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for tumor suppressor research and p53 pathway studies.

[Learn More](#)

Parameter	Value
Product Name	TP53I13 Rabbit pAb
Remarks/Alias	DSCP1, TP53I13
Species	Human
Gene ID	90313 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 180-310 of human TP53I13 (NP_612358.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8NBR0
Protein Weight	42kDa
Shipping	Ice bag