

Anti-Tp53Rk Rabbit Polyclonal Antibody For Wb And Elisa - Q96S44

Item Number: CM0022629



Introduction

TP53RK Rabbit Polyclonal Antibody validated for WB and ELISA; cross-reacts with human, mouse, rat. Targets the EKC/KEOPS complex subunit essential for tRNA modification and p53 regulation. 28 kDa protein. Ideal for IVD and research.

[Learn More](#)

Parameter	Value
Product Name	TP53RK Rabbit pAb
Remarks / Alias	PRPK; BUD32; TPRKB; GAMOS4; Nori-2; Nori-2p; dj101A2; C20orf64; TP53RK
Species	Human
Gene ID (Human)	112858
Gene ID	112858
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-253 of human TP53RK (NP_291028.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96S44
Protein Weight	28kDa
Shipping	Ice bag