

Anti-Tp73 Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - O15350

Item Number: CM0013252



Introduction

Rabbit polyclonal antibody against human TP73 (p73), validated for WB, IF-P, IHC-P, ELISA. Reacts with human and mouse. Targets tumor suppressor p73 involved in apoptosis, cell cycle, and ciliogenesis. Ideal for cancer and neurobiology research.

[Learn More](#)

Parameter	Value
Product Name	TP73 Rabbit pAb
Remarks/Alias	P73; CILD47; TP73
Species	Human
Gene ID (Human)	7161
Gene ID	7161
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 487-636 of human TP73 (NP_005418.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	O15350
Protein Weight	70kDa
Shipping	Ice bag