


Anti-Tpa/Plat Rabbit Monoclonal Antibody For Wb, If/Icc, If-P, Elisa - P00750

Item Number: CM0005623



Introduction

Rabbit monoclonal antibody against human tPA/PLAT, suitable for WB, IF/ICC, IF-P, and ELISA applications. Cross-reacts with mouse and rat. Predicted molecular weight: 63kDa.

[Learn More](#)

Parameter	Value
Product Name	tPA/PLAT Rabbit mAb
Remarks/Alias	TPA; T-PA; tPA/PLAT
Species	Human
GeneID (Human)	5327
GeneID	5327
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 463-562 of human tPA/PLAT (P00750)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00750
Protein Weight	63kDa
Shipping	Ice bag