

Anti-Tpd52 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P55327

Item Number: CM0015706

Introduction

TPD52 rabbit polyclonal antibody (pAb) for Western blot, IHC-P, and ELISA. Targets human tumor protein D52, cross-reacts with mouse and rat. Recombinant immunogen. Ideal for oncology research.

[Learn More](#)

Parameter	Value
Product Name	TPD52 Rabbit pAb
Remarks/Alias	D52; N8L; PC-1; PrLZ; hD52; TPD52
Species	Human
GeneID (Human)	7163
GeneID	7163
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-184 of human TPD52 (NP_005070.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55327
Protein Weight	24kDa
Shipping	Ice bag