

Anti-Tpd52L3 Polyclonal Antibody For Wb, If-P, Elisa - Q96J77

Item Number: CM0015753

Introduction

Rabbit polyclonal antibody targeting human TPD52L3 (D55/hD55/TPD55), validated for Western blot, immunofluorescence, and ELISA. Cross-reacts with mouse and rat. Ideal for testis development studies.

[Learn More](#)

Parameter	Value
Product Name	TPD52L3 Rabbit pAb
Remarks/Alias	D55; hD55; TPD55; NYDSP25; TPD52L3
Species	Human
Gene ID (Human)	89882
Gene ID	89882
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-132 of human TPD52L3 (NP_001001875.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96J77
Protein Weight	16kDa
Shipping	Ice bag