

Anti-Tpd52L3 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - Q96J77

Item Number: CM0017615

Introduction

Rabbit polyclonal antibody against human TPD52L3 (tumor protein D55), validated for WB, IF-P, and ELISA. Cross-reacts with mouse and rat. Suitable for testis-specific protein research.

[Learn More](#)

Parameter	Value
Product Name	TPD52L3 Rabbit pAb
Remarks/Alias	D55, hD55, TPD55, NYDSP25, TPD52L3
Species	Human
Gene ID (Human)	89882
Gene ID	89882
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-132 of human TPD52L3 (NP_001001875.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q96J77
Protein Weight	16kDa
Shipping	Ice bag