


Anti-Tph2 Polyclonal Antibody For Wb And Elisa - Q8Iwu9

Item Number: CM0015525



Introduction

Rabbit polyclonal antibody targeting TPH2 (tryptophan hydroxylase 2), validated for WB and ELISA. Reacts with human, mouse, and rat TPH2. Suitable for serotonin biosynthesis and neurological disorder research.

[Learn More](#)

Specification	Value
Product Name	TPH2 Rabbit pAb
Remarks/Alias	NTPH; ADHD7; TPH2
Species	Human
GeneID (Human)	121278
GeneID	121278
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-70 of human TPH2 (NP_775489.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IWU9
Protein Weight	56kDa
Shipping	Ice bag