

Anti-Tph2 Polyclonal Antibody For Wb, Elisa - Q8Iwu9

Item Number: CM0026952



Introduction

Rabbit polyclonal antibody against TPH2 (neuronal tryptophan hydroxylase) for WB and ELISA, reactive with human, mouse, and rat. Ideal for serotonin biosynthesis and neurological research.

[Learn More](#)

Parameter	Value
Product Name	TPH2 Rabbit pAb
Remarks/Alias	NTPH; ADHD7; TPH2
Species	Human
Gene ID (Human)	121278
Gene ID	121278
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IWU9
Protein Weight	56kDa
Shipping	Ice bag