

Anti-Tpi1 Polyclonal Antibody For Wb, If-P, Elisa - P60174

Item Number: CM0024609



Introduction

Rabbit polyclonal antibody against human triosephosphate isomerase (TPI1), validated for WB, IF-P, and ELISA. Cross-reacts with mouse and rat. Suitable for glycolysis and metabolism research.

[Learn More](#)

Parameter	Value
Product Name	TPI1 Rabbit pAb
Remarks/Alias	TIM; TPI; TPID; HEL-S-49; TPI1
Species	Human
GeneID (Human)	7167
GeneID	7167
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-249 of human TPI1 (NP_000356.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P60174
Protein Weight	27kDa
Shipping	Ice bag