

Anti-Tpi1 Polyclonal Antibody For Wb, Ihc, If/Icc, Elisa - P60174

Item Number: CM0022875



Introduction

High-quality rabbit polyclonal antibody against human TPI1 (Triosephosphate isomerase) validated for WB, IHC, IF/ICC, and ELISA. Detects human, mouse, and rat TPI1.

[Learn More](#)

Parameter	Value
Product Name	TPI1 Rabbit pAb
Remarks/Alias	TIM; TPI; TPID; HEL-S-49; TPI1
Species	Human
Gene ID (Human)	7167
Immunogen	Recombinant Protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P60174
Protein Weight	27kDa
Shipping	Ice bag