

Anti-Tpk1 Polyclonal Antibody For Wb And Elisa - Q9H3S4

Item Number: CM0029159



Introduction

High-quality rabbit polyclonal antibody against human TPK1, validated for Western blot and ELISA. Detects human, mouse, and rat TPK1. Ideal for thiamine metabolism research.

[Learn More](#)

Parameter	Value
Product Name	TPK1 Rabbit pAb
Remarks/Alias	PP20; HTPK1; THMD5; TPK1
Species	Human
Gene ID	27010
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-243 of human TPK1 (NP_071890.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H3S4
Protein Weight	27kDa
Shipping	Ice bag