

Anti-Tpm2 Rabbit Polyclonal Antibody - P07951

Item Number: CM0028367

Introduction

Rabbit polyclonal antibody targeting human TPM2 (Tropomyosin beta chain). Validated for WB, IHC-P, IF/ICC, IP, ELISA applications. Detects TPM2 in human, mouse, and rat tissues. Ideal for muscle biology and cytoskeleton research.

[Learn More](#)

Parameter	Value
Product Name	TPM2 Rabbit pAb
Remarks/Alias	DA1; DA2B; NEM4; TMSB; AMCD1; DA2B4; HEL-S-273; TPM2
Species	Human
GeneID (Human)	7169
GeneID	7169
Immunogen	Synthetic peptide within amino acids 1-100 of human TPM2 (NP_003280.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07951
Protein Weight	33kDa
Shipping	Ice bag