


Anti-Tpmt Rabbit Monoclonal Antibody For Wb And Elisa - P51580

Item Number: CM0002562



Introduction

Rabbit monoclonal antibody targeting human TPMT (thiopurine S-methyltransferase). Designed for WB and ELISA applications, with cross-reactivity to mouse and rat. Enables investigation of thiopurine drug metabolism and pharmacogenomics.

[Learn More](#)

Parameter	Value
Product Name	TPMT Rabbit mAb
Remarks/Alias	TPMTD; TPMT
Species	Human
GeneID (Human)	7172
GeneID	7172
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human TPMT (P51580)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51580
Protein Weight	28kDa
Shipping	Ice bag