


Anti-Tpmt Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P51580

Item Number: CM0016244



Introduction

Rabbit polyclonal antibody against human TPMT (Thiopurine S-methyltransferase), validated in WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for thiopurine drug metabolism research.

[Learn More](#)

Parameter	Value
Product Name	TPMT Rabbit pAb
Remarks/Alias	TPMTD; TPMT
Species	Human
GeneID (Human)	7172
GeneID	7172
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-245 of human TPMT (NP_000358.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51580
Protein Weight	28kDa
Shipping	Ice bag