


Anti-Tpo-R/Cd110/C-Mpl Monoclonal Antibody For Wb, Elisa - P40238

Item Number: CM0003036



Introduction

Rabbit monoclonal antibody targeting human thrombopoietin receptor (TPO-R/CD110/c-Mpl). Validated for WB and ELISA with cross-reactivity to human, mouse, and rat. Essential for research on megakaryocyte differentiation and platelet production.

[Learn More](#)

Parameter	Value
Product Name	TPO-R/CD110/c-Mpl Rabbit mAb
Remarks/Alias	MPLV; TPOR; C-MPL; CD110; THPOR; THCYT2; TPO-R/CD110/c-Mpl
Species	Human
Gene ID (Human)	4352
Gene ID	4352
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 536-635 of human TPO-R/CD110/c-Mpl (P40238).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40238
Protein Weight	71kDa
Shipping	Ice bag