

Anti-Tpp2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P29144

Item Number: CM0025257

Introduction

Rabbit polyclonal antibody against human TPP2 (Tripeptidyl-peptidase 2), validated for WB, IF/ICC, ELISA; cross-reacts with mouse and rat. Ideal for ubiquitin-proteasome pathway and amino acid homeostasis research.

[Learn More](#)

Parameter	Value
Product Name	TPP2 Rabbit pAb
Remarks/Alias	IMD78; TPP-2; TPPII; TPP-II; TPP2
Species	Human
GeneID (Human)	7174
GeneID	7174
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human TPP2 (NP_003282.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P29144
Protein Weight	138kDa
Shipping	Ice bag