

Anti-Tppp3 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Bw30

Item Number: CM0027621



Introduction

Rabbit polyclonal antibody targeting human TPPP3, suitable for WB and ELISA. Cross-reacts with mouse and rat. Ideal for microtubule dynamics, embryo implantation, and endometrial function research. Protein weight 19 kDa.

[Learn More](#)

Parameter	Value
Product Name	TPPP3 Rabbit pAb
Remarks/Alias	p20; CGI-38; TPPP/p20; p25gamma; TPPP3
Species	Human
Gene ID (Human)	51673
Gene ID	51673
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-176 of human TPPP3 (NP_057048.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BW30
Protein Weight	19kDa
Shipping	Ice bag