

Anti-Tppp/P25 Monoclonal Antibody For Wb, Ihc-P, Elisa - O94811

Item Number: CM0008354



Introduction

Rabbit monoclonal antibody targeting human TPPP/p25 (p24, TPPP1). Validated in WB, IHC-P, ELISA; cross-reacts with mouse, rat. For microtubule dynamics, myelination, cell migration, and cell cycle studies. 24kDa.

[Learn More](#)

Parameter	Value
Product Name	TPPP/p25 Rabbit mAb
Remarks/Alias	p24; p25; TPPP1; TPPP/p25; p25alpha
Species	Human
Gene ID (Human)	11076
Gene ID	11076
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 100-200 of human TPPP/p25 (O94811).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O94811
Protein Weight	24kDa
Shipping	Ice bag