

Anti-Tppp/P25 Rabbit Polyclonal Antibody For Wb And Elisa - O94811

Item Number: CM0030459



Introduction

Rabbit polyclonal antibody against TPPP/p25 for Western blot and ELISA. Validated in human, mouse, rat. Suitable for microtubule dynamics and myelination research. - O94811

[Learn More](#)

Parameter	Value
Product Name	TPPP/p25 Rabbit pAb
Remarks/Alias	p24; p25; TPPP1; TPPP/p25; p25alpha
Species	Human
GeneID (Human)	11076
GeneID	11076
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-219 of human TPPP/p25 (NP_008961.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O94811
Protein Weight	24kDa
Shipping	Ice bag