

Anti-Tprg1L Rabbit Polyclonal Antibody For Wb, Elisa - Q5T0D9

Item Number: CM0022360

Introduction

Rabbit polyclonal antibody targeting TPRG1L (mover/FAM79A), a presynaptic protein regulating synaptic release probability. Validated for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for neuroscience research, this antibody detects TPRG1L (MOVER) in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	TPRG1L Rabbit pAb
Remarks/Alias	TPRGL; mover; FAM79A; SVAP30; h-mover; TPRG1L
Species	Human
Gene ID (Human)	127262
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-90 of human TPRG1L (NP_877429.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q5T0D9
Protein Weight	30kDa
Shipping	Ice bag