


Anti-Tpsab1 Rabbit Polyclonal Antibody For Wb And Elisa - Q15661

Item Number: CM0024265



Introduction

Rabbit polyclonal antibody targeting human TPSAB1 (tryptase alpha/beta-1), validated for Western blot and ELISA applications. Cross-reacts with mouse and rat. Suitable for mast cell and innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	TPSAB1 Rabbit pAb
Remarks/Alias	TPS1; TPS2; TPSB1; TPSB2; Tryptase-2; TPSAB1
Species	Human
Gene ID	7177
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 31-275 of human TPSAB1 (NP_003285.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15661
Protein Weight	31kDa
Shipping	Ice bag