


## Anti-Tpst2 Polyclonal Antibody For Wb And Elisa - O60704

Item Number: CM0019769



### Introduction

Rabbit polyclonal antibody against human TPST2, validated for Western blot and ELISA. Detects protein-tyrosine sulfotransferase 2, a Golgi membrane enzyme catalyzing tyrosine O-sulfation. Suitable for human and mouse samples.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | TPST2 Rabbit pAb                                                                                                   |
| Remarks/Alias    | TPST-2; TANGO13B; TPST2                                                                                            |
| Species          | Human                                                                                                              |
| GeneID (human)   | 8459                                                                                                               |
| GeneID           | 8459                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 280-370 of human TPST2 (NP_003586.3) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, ELISA                                                                                                          |
| Cross Reactivity | Human, Mouse                                                                                                       |
| SWISS            | O60704                                                                                                             |
| Protein Weight   | 42kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |