


Anti-Tpx2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9ULW0

Item Number: CM0022663



Introduction

TPX2 rabbit polyclonal antibody, validated for WB, IHC-P, IF/ICC, ELISA; reacts with human and rat. Essential for spindle assembly, mitosis, and cell cycle studies.

[Learn More](#)

Parameter	Value
Product Name	TPX2 Rabbit pAb
Remarks/Alias	DIL2; p100; DIL-2; HCTP4; FLS353; HCA519; REPP86; C20orf1; C20orf2; GD:C20orf1; TPX2
Species	Human
Gene ID	22974
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human TPX2 (NP_036244.2))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9ULW0
Protein Weight	86kDa
Shipping	Ice bag