


Anti-Tra2B Rabbit Monoclonal Antibody For Wb, Ip, If/Icc, Ihc-P, Elisa - P62995

Item Number: CM0012397



Introduction

This anti-TRA2B rabbit monoclonal antibody is offered for WB, IP, IF/ICC, IHC-P, and ELISA, cross-reacting with human, mouse, and rat. Ideal for studies on pre-mRNA splicing regulation and neurodegenerative diseases.

[Learn More](#)

Parameter	Value
Product Name	TRA2B Rabbit mAb
Remarks/Alias	SFRS10; SRF510; TRAN2B; PPP1R156; TRA2-BETA; Htra2-beta
Species	Human
GeneID (Human)	6434
GeneID	6434
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human TRA2B (NP_004584.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62995
Protein Weight	4 kDa/22 kDa/34 kDa
Shipping	Ice bag