


Anti-Tra2B Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P62995

Item Number: CM0026731



Introduction

Anti-TRA2B rabbit polyclonal antibody validated for WB, IF/ICC, ELISA. Recognizes human, mouse, rat TRA2B. Ideal for splicing regulation and neurodegeneration studies. - P62995

[Learn More](#)

Parameter	Value
Product Name	TRA2B Rabbit pAb
Remarks/Alias	SFRS10; SRF510; TRAN2B; PPP1R156; TRA2-BETA; Htra2-beta; TRA2B
Species	Human
GeneID	6434
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human TRA2B (NP_004584.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62995
Protein Weight	34kDa
Shipping	Ice bag