

Anti-Trac Rabbit Polyclonal Antibody For Wb And Elisa - P01848

Item Number: CM0019877



Introduction

Rabbit polyclonal antibody targeting human T cell receptor alpha chain constant region (TRAC). Validated for Western blot and ELISA, with cross-reactivity to mouse. Suitable for T cell immunology and immunoassay development.

[Learn More](#)

Parameter	Value
Product Name	TRAC Rabbit pAb
Remarks/Alias	IMD7; TCRA; TRA@; TRAC
Species	Human
GeneID (Human)	6955
GeneID	6955
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-176 of human TRAC (AAH20840.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P01848
Protein Weight	16kDa
Shipping	Ice bag