

Anti-Tradd Rabbit Polyclonal Antibody For Wb & Elisa - Q15628

Item Number: CM0014429



Introduction

High-quality anti-TRADD rabbit polyclonal antibody validated for Western blot and ELISA applications. Detects human and mouse TRADD, a key adaptor in TNF-induced apoptosis and NF- κ B activation.

[Learn More](#)

Parameter	Value
Product Name	TRADD Rabbit pAb
Remarks/Alias	9130005N23Rik; TRADD
Species	Human
GeneID (human)	8717
GeneID	8717
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 199-312 of human TRADD (NP_003780.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15628
Protein Weight	34kDa
Shipping	Ice bag