


Anti-Traf2 Monoclonal Antibody (Kd Validated) For Wb, Ip, Elisa - Q12933

Item Number: CM0002653



Introduction

Rabbit monoclonal [KD Validated] TRAF2 antibody for WB, IP, ELISA. Reacts with human, mouse, rat. Targets E3 ubiquitin ligase involved in NF- κ B and JNK signaling, critical for apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] TRAF2 Rabbit mAb
Remarks/Alias	TRAP; TRAP3; RNF117; MGC:45012; [KD Validated] TRAF2
Species	Human
GeneID (Human)	7186
GeneID	7186
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12933
Protein Weight	53kDa/54kDa/55kDa/61kDa
Shipping	Ice bag