


Anti-Traf2 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q12933

Item Number: CM0014192



Introduction

Rabbit polyclonal antibody against human TRAF2 (SWISS: Q12933). Validated for WB, IHC-P, ELISA; cross-reacts with human, mouse, rat. TRAF2 is an E3 ubiquitin ligase regulating NF- κ B, JNK, and apoptosis pathways. Ideal for signal transduction research.

[Learn More](#)

Parameter	Value
Product Name	TRAF2 Rabbit pAb
Remarks/Alias	TRAP; TRAP3; RNF117; MGC:45012; TRAF2
Species	Human
Gene ID (Human)	7186
Gene ID	7186
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human TRAF2 (NP_066961.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12933
Protein Weight	56kDa
Shipping	Ice bag