


Anti-Traf2 Rabbit Polyclonal Antibody For Wb And Elisa - Q12933

Item Number: CM0020289



Introduction

TRAF2 rabbit polyclonal antibody validated for WB and ELISA, reactive with human, mouse, and rat. Targets TNF receptor-associated factor 2, a key regulator of NF- κ B and JNK signaling involved in cell survival and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	TRAF2 Rabbit pAb
Remarks/Alias	TRAP; TRAP3; RNF117; MGC:45012; TRAF2
Species	Human
GeneID (Human)	7186
GeneID	7186
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 190-350 of human TRAF2 (NP_066961.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12933
Protein Weight	56 kDa
Shipping	Ice bag